



## Multifunction LED light system for Face and Body treatment

- Anti-Ageing;
- Bio stimulation
- Toning
- Drainage
- Activation of active principles



## Sistema multifunzionale a Luce LED per trattamenti Viso e Corpo:

- Anti-Ageing;
- Bio stimolazione
- Tonificazione
- Drenaggio
- Veicolazione di Principi Attivi




ER  
E  
S  
A  
U  
Q



**QUASER** is a new generation light emitting diodes machine for body and face treatment. It employs the principles of photo dynamic therapy for anti-ageing treatment , bio stimulation, drainage and activation of active principles.

**QUASER** è un apparecchiatura a diodi led di nuova generazione per il trattamento del viso e del corpo. Si avvale dei principi della terapia fotodinamica per trattamenti anti-ageing , biostimolazione, drenaggio ed attivazione di principi attivi

CARATTERISTICHE TECNICHE		TECHNICAL SPECIFICATIONS
ALIMENTAZIONE	230V—50/60 Hz	MAIN SUPPLY
POTENZA ASSORBITA	120 VA	ABSORBED POWER
PESO (kg)	7 kg	WEIGHT (kg)
DIMENSIONI	360x360x320 mm	DIMENSION
RGB	72 RGB sources 0,5W each Rosso / Red: 625nm Verde / Green: 525nm Blu / Blue: 465nm	RGB
FLUENZA TOTALE	700 Lu W	TOTAL FLUENCY
RIPETIZIONE IMPULSO	2-10Hz	PULSE REPETITION
AREA ZONA VISO	600cm <sup>2</sup>	TREATMENT AREA(FACE)
AREA ZONA CORPO	600 cm <sup>2</sup>	TREATMENT AREA(BODY)
CONTROLLO	Microprocessore—Microprocessor	ELECTRONC CONTROL



**TecnoLaser**

Via Venezia 173 - 43122 Parma - ITALY  
Tel. +39 0521 271102  
Fax +39 0521 272846  
e-mail: [tecnolaser@tin.it](mailto:tecnolaser@tin.it)

[www.tecnolaser.net](http://www.tecnolaser.net) / [www.tecnolaser.de/](http://www.tecnolaser.de/)  
[www.technolaser.eu/](http://www.technolaser.eu/) / [www.tecnolasereuropa.info/](http://www.tecnolasereuropa.info/)  
[www.tecnogamma.info/](http://www.tecnogamma.info/) / [www.applegreen.net](http://www.applegreen.net)  
[www.applegreen.it/](http://www.applegreen.it/) / [www.applegreen.eu](http://www.applegreen.eu)