





Multifunction LED light system for Face and Body treatment

- Anti-Ageing;
- Bio stimulation
- Toning
- Drainage
- Activation of active principles



Sistema multifunzionale a Luce LED per trattamenti Viso e Corpo:

- Anti-Ageing;
- Bio stimolazione
- Tonificazione
- Drenaggio
- Veicolazione di Principi Attivi




ER
E
S
A
U
Q



QUASER is a new generation light emitting diodes machine for body and face treatment. It employs the principles of photo dynamic therapy for anti-ageing treatment , bio stimulation, drainage and activation of active principles.

QUASER è un apparecchiatura a diodi led di nuova generazione per il trattamento del viso e del corpo. Si avvale dei principi della terapia fotodinamica per trattamenti anti-ageing , biostimolazione, drenaggio ed attivazione di principi attivi

CARATTERISTICHE TECNICHE		TECHNICAL SPECIFICATIONS
ALIMENTAZIONE	230V—50/60 Hz	MAIN SUPPLY
POTENZA ASSORBITA	120 VA	ABSORBED POWER
PESO (kg)	7 kg	WEIGHT (kg)
DIMENSIONI	360x360x320 mm	DIMENSION
RGB	72 RGB sources 0,5W each Rosso / Red: 625nm Verde / Green: 525nm Blu / Blue: 465nm	RGB
FLUENZA TOTALE	700 Lu W	TOTAL FLUENCY
RIPETIZIONE IMPULSO	2-10Hz	PULSE REPETITION
AREA ZONA VISO	600cm ²	TREATMENT AREA(FACE)
AREA ZONA CORPO	600 cm ²	TREATMENT AREA(BODY)
CONTROLLO	Microprocessore—Microprocessor	ELECTRONC CONTROL



TecnoLaser

Via Venezia 173 - 43122 Parma - ITALY
Tel. +39 0521 271102
Fax +39 0521 272846
e-mail: tecnolaser@tin.it

www.tecnolaser.net / www.tecnolaser.de/
www.technolaser.eu/ / www.tecnolasereuropa.info/
www.tecnogamma.info/ / www.applegreen.net
www.applegreen.it/ / www.applegreen.eu